




ROI Manager

Region Of Interest (ROI)




Secondary/Right Click
for more options



Double Click to
set/change properties

Straight, segmented or freehand lines, or arrows (right c


ROI Manager



Analyze > Tools > ROI Manager...

† *


The ROI Manager is a tool for working with **multiple** selections (ROIs).



*(cmd) + t can be used to both to open the ROI Manager and/or add a new ROI to the Manager.

**shift + e can be used to draw the last ROI.

ROI Manager



Add new ROI (t).




After modifying a ROI, you can use Update to save the changes.

Delete selected ROI. If none is selected, delete all.



Rename selected ROI.

Deselect one or more selected ROI.

ROI Manager





Change the properties of the selected ROI(s)



“Analyze” menu

“Edit” menu



ROI Manager



Open a saved ROI.zip file (drag & drop works as well)


Save ROI(s) as a .zip file

ROI Manager



Open a saved ROI.zip file (drag & drop works as well)
Save ROI(s) as a .zip file

Set ROI(s) label options



*Labels option can also be found under “Image > Overlay”


*“Use names as labels” can be also found under “More>Options...”




Live Demo of ROI manager

Segmentation And Measurements




Binary mask





Process and segment





Export and Measure




Binary mask




Process and segment



Export and Measure




Do you see an issue here?




in **Fiji**: **Process** > **Binary** > **Watershed**

Watershed is a useful algorithm to try to **separate touching objects**.


Image




Binary Mask







Watershed





Binary mask




Distance transform






Binary mask




Instance segmentation




in **Fiji**: **Analyze > Analyze Particles...**




in **Fiji**:
Analyze
> **Analyze Particles...**




in **Fiji**:
Analyze
> **Analyze Particles...**




in **Fiji**:
Analyze
> **Analyze Particles...**



in **Fiji**:
Analyze
> **Analyze Particles...**




in **Fiji**:
Analyze
> **Analyze Particles...**




Analyze particles: Select for morphology




in Fiji:
Analyze
> **Analyze Particles...**





Size




Circularity




Fill




in **Fiji**: **Analyze > Analyze Particles...**




Size range of the **particles** that you want to detect.






the **size range** of the **particles** that you want to detect.




in **Fiji**: **Analyze > Analyze Particles...**



Circularity of the **particles** that you want to detect.





$$4\pi \times \frac{\text{area}}{\text{perimeter}^2} = 1$$



Analyze particles: Circularity






$$4\pi \times \frac{\text{area}}{\text{perimeter}^2} = 1$$




$$4\pi \times \frac{\text{area}}{\text{perimeter}^2} < 1$$



$$4\pi \times \frac{\text{area}}{\text{perimeter}^2} = 1$$




in **Fiji**: **Analyze > Analyze Particles...**




“**Exclude on edges**” excludes objects that are touching the borders of the image.






Binary mask



Process and segment




Export and Measure







- Save the segmented image
- Add to ROI manager
- Generate and export measurements

Export and Measure




Results			
	Area	Mean	X
0016	285	204.29474	197.
0028	278	174.84892	219.
0027	231	188.46753	45.0
0034	501	189.14172	174.
0041	660	171.69697	73.8
0048	228	195.89474	233.
0055	448	209.03571	138.

in **Fiji**: **Analyze > Analyze Particles...**




“**Count Masks**” encodes object identity as gray value.





in Fiji: **Analyze > Analyze Particles...**



“Add to Manager” stores all found objects in the ROI Manager.




in **Fiji**: **Analyze > Analyze Particles...**



- Save the segmented image
- Add to ROI manager

- Generate and export measurements

Export and Measure




The image shows a grayscale microscopy image of cells. The cells are segmented and labeled with numbers from 1 to 37. A table of measurements is overlaid on the image, showing the results for each cell. The table has four columns: ID, Area, Mean, and X. The data is as follows:

Results			
	Area	Mean	X
0016	285	204.29474	197.
0028	278	174.84892	219.
0027	231	188.46753	45.0
0034	501	189.14172	174.
0041	660	171.69697	73.8
0048	228	195.89474	233.
0055	448	209.03571	138.

- Save the segmented image
- Add to ROI manager

- Generate and export measurements
 - Select what to measure
 - Measure

Export and Measure





The image shows a grayscale microscopy image of cells. The cells are segmented and labeled with numbers from 1 to 39. A table of measurements is overlaid on the image, showing the results for each cell. The table has four columns: ID, Area, Mean, and X. The data is as follows:

Results			
	Area	Mean	X
0016	285	204.29474	197.
0028	278	174.84892	219.
0027	231	188.46753	45.0
0034	501	189.14172	174.
0041	660	171.69697	73.8
0048	228	195.89474	233.
0055	448	209.03571	138.

Analyze > Set Measurements...

Specifies which measurements have to be performed (e.g. area, mean grey value, max and min grey values, ...)








Analyze > Measure

(cmd) + m



Measures the parameters chosen under “Analyze > Set Measurements...” in relation to the selected ROI.

Results are displayed in a Result Table (which can be saved as .csv, .xlsx, ...)







Analyze > Plot Profile

(cmd) + k



*Plots can be saved as .csv file ("Save Data...") and also as images (e.g. "File > Save AS > PNG")





1. Select the “ORIGINAL” image.
2. To measure all the ROIs stored in the “ROI Manager”, **no ROI has to be selected** in the “ROI Manager”. Use the “Deselect” button.
3. Click on the “**Measure**” button.
4. A “**Results**” table will appear.

Segmentation with pixel based classifier—exercises



Continue with the “Analyze Particle” step from ONE of the workflow exercises below.

4.1 DAPI segmentation with thresholding

OR

4.3 - DAPI segmentation with Labkit